F&P Evora Full | Sizing guide


Printing instructions

- Print the sizing guide
 - Measure your printout > Place your ruler on the scale below and ensure they measure the same.
- > Print on U.S. Letter- or A4-sized paper.
- > Ensure you select "Actual size" and not "Fit" in the printer settings.
- > For best results, print in color.

6 7 8 9 10 11 12 13 14

Sizing guide instructions

Carefully cut along the dotted lines.



Hold the guide so that it presses against the underside of your nose.



Choose the size that best fits your nose.

